



## A Comprehensive Solution for HPV Genotyping

Distinguish 14 types of HPV in a single tube




14 high-risk HPV genotypes

14 low-risk HPV genotypes

**HPV** is the most common sexually transmitted infections(STI). There were about 43 million HPV infections in 2018, many of them were found in their late teens and early 20s. There are many different types of HPV. Some types can cause health problems, including genital warts and cancers.

Overall HPV prevalence and percentage of each genotype in different regions



## HPV Product Advantages



Single-tube for **14** high-risk or **14** low-risk HPV genotyping



**94** samples in **2.5** hours



Patented technology **MMCA**<sup>®</sup>



Automated **Result** interpretation

## Reagent Information

| Items          | Information  |
|----------------|--|
| Types          | 14 High-Risk HPV: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68<br>14 Low-Risk HPV: 6, 11, 26, 40, 42, 43, 44, 53, 54, 61, 69, 70, 72, 81 |
| Reagent type   | Pre-loaded freeze-dry reagent  |
| Sample size    | 25 µL  |
| Time to result | 2.5 hours  |
| LoD            | 500 copies/reaction  |
| Specificity    | 98.63%   |
| Reference gene | GAPDH  |
| Result         | Automated interpretation   |



# HPV Genotyping Detection Workflow

## Option 1 → MeltPro® HPV Genotyping Detection Workflow

### Automated High Throughput Extraction Solutions (96 samples)

Lab-Aid 896 Nucleic Acid Extraction System



process up to **96** samples in a single run, in batches of **12**


### Automated Medium Throughput Extraction (24 samples)

Lab-Aid 824s Nucleic Acid Extraction System




process **1-24** samples in a single run

### Multicolor Melting Curve Analysis (MMCA®)



### Sample Report (e.g. HPV 16 positive)



## Option 2 → Integrated Detection Platform

01

Sample Loading




02

Sanity 2.0 System





## ● Screening for 14 types of high-risk HPV



## ● Verification of experimental results

was originally found to be 56%, but, when adjusted, this accuracy rose to 80% (Table 2).

**DISCUSSION**

Near-patient hrHPV testing in LMICs has the potential to improve outcomes of cervical cancer screening programs by reducing screening costs with use of high-throughput PCR devices, improving care efficiency, and improving retention of at-risk participants. The ideal LMIC test should be sensitive, specific, user friendly, and robust; require minimal equipment; and be rapid and affordable.<sup>21</sup> With those attributes in mind, hrHPV MMCA near-patient testing was implemented in two separate field locations, each providing representative challenges typical of LMICs.

The hrHPV MMCA assay was highly sensitive and specific at both field sites and accuracy of HPV genotyping enabled additional participant risk stratification, potentially reducing the impact of screening on a taxed health-care system. Since the link of cervical cancer to HPV infection was made, approximately 70% of cervical cancer cases have been linked to infection with HPV16 or

box tests to avoid multiple procedural steps. The hrHPV MMCA achieves this by providing lyophilized hrHPV master mix prepackaged into PCR strip tubes. Temperatures at field sites in Honduras reached > 32.2°C, with high levels of humidity and no access to refrigeration for reagent storage. After 4 days of exposure to elevated temperatures, no differences in test sensitivity or specificity were observed. In addition, minimally trained personnel needed only to pipet crude lysate directly into the lyophilized reagents and place the reaction tubes into the quantitative PCR (qPCR) thermocycler for successful analysis. As noted, invalid rates were highest during the testing at the first field site but dropped in the more infrastructure-challenged

**Table 2. Genotyping Accuracy of Multicolor Melt Curve Analysis Using Field Protocol**

| HPV Type    | Accuracy (%) | Adjusted Accuracy (%) |
|-------------|--------------|-----------------------|
| 16          | 100          | 100                   |
| 18          | 100          | 100                   |
| 59          | 78           | 78                    |
| 68          | 88           | 100                   |
| 58          | 91           | 100                   |
| 31          | 100          | 100                   |
| 39          | 100          | 100                   |
| Confections | 56           | 80                    |

NOTE. Out-of-the-box multicolor melt curve analysis was used to

MeltPro® HPV offers high sensitivity and specificity, which reduces the risk of progression to cervical precancer / cancer due to undefined genotypes.

Experimental results shows the high stability of MeltPro® HPV.

Turner S A , Deharvengt S J , Lyons K D , et al. Implementation of Multicolor Melt Curve Analysis for High-Risk Human Papilloma Virus Detection in Low- and Middle-Income Countries: A Pilot Study for Expanded Cervical Cancer Screening in Honduras[J]. Journal of Global Oncology, 2018(4):1-8.

## Products Information

| Cat.No. | Product name                                   | Specification / Model |
|---------|--|-----------------------|
| 604106  | Lab-Aid 824 HPV DNA Extraction Kit             | 48 Tests/Kit          |
| 610104  | HPV DNA Extraction Kit (Lab-Aid 896)           | 96 Tests/Kit          |
| 608112  | HPV DNA Extraction Kit (Sanity 2.0)            | 24 Tests/Kit          |
| 801170  | High Risk HPV Genotyping Kit (Sanity 2.0)      | 24 Tests/Kit          |
| 801011  | MeltPro® High Risk HPV Genotyping Assay (MMCA) | 48 Tests/Kit          |
| 801013  | MeltPro® Low Risk HPV Genotyping Assay         | 48 Tests/Kit          |
| 505101  | Sanity 2.0 System                              | Sanity® 2.0           |
| 502602  | Real-Time PCR System (SLAN-48P)                | SLAN® 48P             |
| 502601  | Real-Time PCR System (SLAN-96P)                | SLAN® 96P             |
| 502603  | Real-Time PCR System (SLAN-96S)                | SLAN® 96S             |
| 501101  | Lab-Aid 824s Nucleic Acid Extraction System    | Lab-Aid® 824s         |
| 501110  | Nucleic Acid Extraction System (Lab-Aid 896)   | Lab-Aid® 896          |

For more information, please visit <http://www.zeesandx.com>

### Xiamen Zeesan Biotech Co., Ltd.

Building 1#, 4F and A area 3F of Building 11#, No. 3701 North Xiang'an Road, (Xiang An) Industrial Area, Torch High-Tech Zone, Xiamen City, Fujian Province, 361101, P.R. China

Phone: +86-592-7615091

E-mail: info@zsandx.com